

## Swim tunnel specifications

Specifications	VOLUME [L]	TEST SECTION [cm]	FISH SIZE [g]	SPEED [cm/s]	FOOT PRINT [cm]	MOTOR SIZE [kW]	POWER	FLUSH PUMP	FLUSH FLOW [L/min]
SW10000	0,17	ID2,64 x L10	1-4	0,7-50	75 x 25	0,046	230 V / 50 Hz / 0,5 A	PU10050	5
SW10020	0,17	ID2,64 x L10	1-4	0,7-50	75 x 25	0,046	110 V / 60 Hz / 1,0 A	PU10060	5
SW10030	1,5	ID5,5 x L20	4-12	3-100	136 x 25	0,046	230 V / 50 Hz / 0,5 A	PU10002	5
SW10040	1,5	ID5,5 x L20	4-12	3-100	136 x 25	0,046	110 V / 60 Hz / 1,0 A	PU10060	5
SW10050	5	30 x 7,5 x 7,5	20-80	3-110	117 x 40	0,37	230 V / 50 Hz / 3,2 A	PU10100	10
SW10060	5	30 x 7,5 x 7,5	20-80	3-110	117 x 40	0,37	110 V / 60 Hz / 6,4 A	PU10110	10
SW10100	10	40 x 10 x 10	50-150	3-110	130 x 45	0,37	230 V / 50 Hz / 3,2 A	PU10150	20
SW10110	10	40 x 10 x 10	50-150	3-110	130 x 45	0,37	110 V / 60 Hz / 6,4 A	PU10160	20
SW10150	30	55 x 14 x 14	175-500	3-110	147 x 53	0,37	230 V / 50 Hz / 3,2 A	PU10200	40
SW10160	30	55 x 14 x 14	175-500	3-110	147 x 53	0,37	110 V / 60 Hz / 6,4 A	PU10210	40
SW10200	90	70 x 20 x 20	450-1500	5-150	188 x 71	0,55	230 V / 50 Hz / 4,8 A	PU10250	57
SW10210	90	70 x 20 x 20	450-1500	5-150	188 x 71	0,55	110 V / 60 Hz / 9,6 A	PU10260	57
SW10250	185	87,5 x 25 x 25	750-5000	10-225	227 x 91	4,0	3X360-440 V / 50 Hz / 8,8 A	not included	-
SW10275	185	87,5 x 25 x 25	750-5000	10-225	227 x 91	4,0	3X400-480 V / 60 Hz / 8,0 A	not included	-
SW10300	850	135 x 45 x 45	5000-15000	15-200	132 x 347	11,0	3X360-440 V / 50 Hz / 30,0 A	not included	-
SW10325	850	135 x 45 x 45	5000-15000	15-200	132 x 347	11,0	3X400-480 V / 60 Hz / 30,0 A	not included	-



[www.loligosystems.com](http://www.loligosystems.com)